



GPX7 Antibody

Product Code	CSB-PA009872GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SL4
Immunogen	Human GPX7
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	glutathione peroxidase 7;GPX7;CL683;FLJ14777;GPX6;NPGPx ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPX7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.