



GRIPAP1 Antibody

Product Code	CSB-PA009923GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4V328
Immunogen	Human GRIPAP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	GRIP1 associated protein 1;GRIPAP1;DKFZp434P0630;GRASP-1;KIAA1167;MGC126593;MGC126595 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRIPAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.