



GSPT1 Antibody

Product Code	CSB-PA009966GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15170
Immunogen	Human GSPT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	G1 to S phase transition 1;GSPT1;551G9.2;ETF3A;FLJ38048;FLJ39067;GST1;eRF3a ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GSPT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.