



GTF2H1 Antibody

Product Code	CSB-PA010007GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32780
Immunogen	Human GTF2H1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	general transcription factor IIH, polypeptide 1, 62kDa; GTF2H1; BTF2; TFB1; TFIIH ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTF2H1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.