



GTF3C5 Antibody

Product Code	CSB-PA010029GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Q8
Immunogen	Human GTF3C5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	general transcription factor IIIC, polypeptide 5, 63kDa;GTF3C5;FLJ20857;TFIIIC63;TFIIICepsilon;TFIIC2-63 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTF3C5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.