



GUCY1A1 Antibody

Product Code	CSB-PA010050GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02108
Immunogen	Human GUCY1A3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	guanylate cyclase 1, soluble, alpha 3;GUCY1A3;GC-SA3;GUC1A3;GUCA3;GUCSA3;GUCY1A1 ;
Immunogen Species	Homo sapiens (Human)
Target Names	GUCY1A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.