



HADH Antibody

Product Code	CSB-PA010117GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16836
Immunogen	Human HADH
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	hydroxyacyl-Coenzyme A dehydrogenase;HADH;HAD;HADH1;HADHSC;HHF4;M/SCHAD;MGC8392;SCHAD ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HADH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.