



HAVCR2 Antibody

Product Code	CSB-PA010145GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDQ0
Immunogen	Human HAVCR2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	hepatitis A virus cellular receptor 2; HAVCR2; FLJ14428; KIM-3; TIM3; TIMD3; Tim-3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HAVCR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.