



HDHD3 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA010256GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BSH5 |
| Immunogen | Human HDHD3 |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | haloacid dehalogenase-like hydrolase domain containing 3; HDHD3; 2810435D12Rik; C9orf158; MGC12904 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HDHD3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |