



Recombinant Pseudomonas mendocina Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033PZM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4XPX3
Product Type	Recombinant Protein
Immunogen Species	Pseudomonas mendocina (strain ymp)
Purity	>85% (SDS-PAGE)
Sequence	MKQRLFILSQ YLLPHHLLSR LIGCVAECRI GWLKNRLIAW FAKQYRVDMS EAQVEDLSAY EHFNAFFTRA LKAGARPLDE TPGAILCPAD GAISQLGVIE HGRLFQAKGH SFSATELLGG DAERAAPFMG GQFATVYLSP KDYHRVHMPL AGTLKEMVYV PGRLFSVNQT TAENVPELFA RNERVVCLFD TERGPMVAVL VGAMIVASVE TVWAGLVTPP KRELKSTRYD AESRGIPIELA KGAE LGRFKL G
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.