



HIGD1A Antibody

Product Code	CSB-PA010354GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y241
Immunogen	Human HIGD1A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	HIG1 hypoxia inducible domain family, member 1A;HIGD1A;DKFZp564K247;HIG1 ;
Immunogen Species	Homo sapiens (Human)
Target Names	HIGD1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.