



HIP1R Antibody

Product Code	CSB-PA010366GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75146
Immunogen	Human HIP1R
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	huntingtin interacting protein 1 related;HIP1R;FLJ14000;FLJ27022;HIP12;HIP3;ILWEQ;KIAA0655;MGC47513 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIP1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.