



HMGB1 Antibody

Product Code	CSB-PA010553GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09429
Immunogen	Human HMGB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	high-mobility group box 1;HMGB1;DKFZp686A04236;HMG1;HMG3;SBP-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HMGB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.