



HMGN4 Antibody

Product Code	CSB-PA010576GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00479
Immunogen	Human HMGN4
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	high mobility group nucleosomal binding domain 4;HMGN4;HMG17L3;MGC5145;NHC ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HMGN4
Target Details	This protein, a member of the HMGN protein family, is thought to reduce the compactness of the chromatin fiber in nucleosomes, thereby enhancing transcription from chromatin templates. Transcript variants utilizing alternative polyadenylation signals exist for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.