



HPN Antibody

Product Code	CSB-PA010704GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05981
Immunogen	Human Hpn
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Purified Protein A
Isotype	IgG
Alias	hepsin;Hpn;AA408042;AI788568 ;
Product Type	Purified Rabbit Anti Mus PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HPN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.