



HPSE Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA010716GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y251 |
| Immunogen | Human HPSE |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,IHC |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | heparanase;HPSE;HPA;HPR1;HPSE1;HSE1 ; |
| Product Type | Purified Rabbit Anti Human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HPSE |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |