



HRASLS5 Antibody

Product Code	CSB-PA010729GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KN8
Immunogen	Human HRASLS5
Species Reactivity	Human
Tested Applications	ELISA
Alias	;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HRASLS5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.