



HSD17B11 Antibody

Product Code	CSB-PA010768GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBQ5
Immunogen	Human HSD17B11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	hydroxysteroid (17-beta) dehydrogenase 11;HSD17B11;17-BETA-HSD11;17-BETA-HSDXI;DHRS8;PAN1B;RETSR2;SDR16C2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSD17B11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.