



HSP90AB1 Antibody

Product Code	CSB-PA010808GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08238
Immunogen	Human HSP90AB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	heat shock protein 90kDa alpha (cytosolic), class B member 1;HSP90AB1;D6S182;FLJ26984;HSP90-BETA;HSP90B;HSPC2;HSPCB ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSP90AB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.