



HSPA1A Antibody

Product Code	CSB-PA01082GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DMV8
Immunogen	Human HSP70
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	heat shock 70kDa protein 4;HSPA4;APG-2;HS24/P52;MGC131852;RY;hsp70;hsp70RY ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSPA1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.