



HSPH1 Antibody

Product Code	CSB-PA010869GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92598
Immunogen	Human HSPH1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	heat shock 105kDa/110kDa protein 1;HSPH1;DKFZp686M05240;HSP105;HSP105A;HSP105B;KIAA0201;NY-CO-25 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSPH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.