



HYLS1 Antibody

Product Code	CSB-PA010927GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96M11
Immunogen	Human HYLS1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	hydrolethalus syndrome 1; HYLS1; FLJ32915; HLS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HYLS1
Target Details	This gene encodes a protein localized to the cytoplasm. Mutations in this gene are associated with hydrolethalus syndrome. Multiple alternatively spliced variants, encoding the same protein, have been identified.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.