



IFIT3 Antibody

Product Code	CSB-PA011022GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14879
Immunogen	Human IFIT3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	interferon-induced protein with tetratricopeptide repeats 3; IFIT3; CIG-49; GARG-49; IFI60; IFIT4; IRG2; ISG60; RIG-G ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFIT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.