



IFITM1 Antibody

Product Code	CSB-PA011024GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13164
Immunogen	Human IFITM1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	interferon induced transmembrane protein 1 (9-27);IFITM1;9-27;CD225;IFI17;LEU13 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFITM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.