



IFT140 Antibody

Product Code	CSB-PA011072GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RY7
Immunogen	Human IFT140
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	intraflagellar transport 140 homolog (Chlamydomonas);IFT140;DKFZp564L232;FLJ10306;FLJ30571;KIAA0590;WDTC2;c305C8.4;c380F5.1;gs114 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFT140
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.