



IGF2BP2 Antibody

Product Code	CSB-PA011091GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6M1
Immunogen	Human IGF2BP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	insulin-like growth factor 2 mRNA binding protein 2;IGF2BP2;IMP-2;IMP2;VICKZ2;p62 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IGF2BP2
Target Details	This gene encodes a member of the IGF-II mRNA-binding protein (IMP) family. This protein contains several four KH domains and two RRM domains. It functions by binding to the 5' UTR of the insulin-like growth factor 2 (IGF2) mRNA and regulating IGF2 translation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.