



IGHA1 Antibody

Product Code	CSB-PA011111GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01876
Immunogen	Human IGHA1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	immunoglobulin heavy constant alpha 1; IGHA1; FLJ14473; FLJ35065; FLJ35500; FLJ36402; FLJ39698; FLJ40001; FLJ41548; FLJ41552; FLJ41789; FLJ43248; FLJ43594; FLJ44293; FLJ46028; FLJ46621; FLJ46724; FLJ46811; FLJ46824; FLJ90170; IgA1; MGC102857 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IGHA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.