



IGHA1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA011111GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P01876 |
| Immunogen | Human IGHA1 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | immunoglobulin heavy constant alpha 1; IGHA1; FLJ14473; FLJ35065; FLJ35500; FLJ36402; FLJ39698; FLJ40001; FLJ41548; FLJ41552; FLJ41789; FLJ43248; FLJ43594; FLJ44293; FLJ46028; FLJ46621; FLJ46724; FLJ46811; FLJ46824; FLJ90170; IgA1; MGC102857 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IGHA1 |