



Recombinant Human Thioredoxin-like protein 4A (TXNL4A)

Product Code	CSB-YP025386HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P83876
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSYMLPHLHN GWQVDQAILS EEDRVVVIRF GHDWDPTCMK MDEVLYSIAE KVKNFAVIYL VDITEVPDFN KMYELYDPCT VMFFFRNKHI MIDLGTGNNN KINWAMEDKQ EMVDIETVY RGARKGRGLV VSPKDYSTKY RY
Source	Yeast
Target Names	TXNL4A
Protein Names	Recommended name: Thioredoxin-like protein 4A Alternative name(s): DIM1 protein homolog Spliceosomal U5 snRNP-specific 15 kDa protein Thioredoxin-like U5 snRNP protein U5-15kD
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.