



IL20RB Antibody

Product Code	CSB-PA011632GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXL0
Immunogen	Human IL20RB
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	interleukin 20 receptor beta; IL20RB; DIRS1; FNDC6; IL-20R2; MGC34923 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL20RB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.