



IL22 Antibody

Product Code	CSB-PA011635GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZX6
Immunogen	Human IL22
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	interleukin 22;IL22;IL-21;IL-22;IL-D110;IL-TIF;IL21;ILTIF;MGC79382;MGC79384;TIFIL-23;TIFa;zcyto18 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL22
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.