



INTS3 Antibody

Product Code	CSB-PA011761GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68E01
Immunogen	Human INTS3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	integrator complex subunit 3;INTS3;C1orf193;C1orf60;DKFZp686E1950;DKFZp686O20115;DKFZp78111253;FLJ21919;FLJ44849;INT3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	INTS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.