



IPPK Antibody

Product Code	CSB-PA011783GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8X2
Immunogen	Human IPPK
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	inositol 1,3,4,5,6-pentakisphosphate 2-kinase;IPPK;C9orf12;FLJ13163;INSP5K2;IP5K;IPK1;KIAA0699;bA476B13.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IPPK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.