



IRF2BP1 Antibody

Product Code	CSB-PA011816GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU81
Immunogen	Human IRF2BP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	interferon regulatory factor 2 binding protein 1;IRF2BP1;DKFZp434M154 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IRF2BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.