



ISG15 Antibody

Product Code	CSB-PA011843GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05161
Immunogen	Human ISG15
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ISG15 ubiquitin-like modifier;ISG15;G1P2;IFI15;UCRP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ISG15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.