



ITGB1BP1 Antibody

Product Code	CSB-PA011881GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14713
Immunogen	Human ITGB1BP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	integrin beta 1 binding protein 1;ITGB1BP1;DKFZp686K08158;ICAP-1A;ICAP-1B;ICAP-1alpha;ICAP1;ICAP1A;ICAP1B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITGB1BP1
Target Details	The cytoplasmic domains of integrins are essential for cell adhesion. This protein binds to the beta1 integrin cytoplasmic domain. The interaction between this protein and beta1 integrin is highly specific. Two isoforms of this protein are derived from alternatively spliced transcripts. The shorter form of this protein does not interact with the beta1 integrin cytoplasmic domain. The longer form is a phosphoprotein and the extent of its phosphorylation is regulated by the cell-matrix interaction, suggesting an important role of this protein during integrin-dependent cell adhesion.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.