



ITPKC Antibody

Product Code	CSB-PA011911GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DU7
Immunogen	Human ITPKC
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	inositol 1,4,5-trisphosphate 3-kinase C;ITPKC;IP3KC ;
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITPKC