



IVNS1ABP Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA011923GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y6Y0 |
| Immunogen | Human IVNS1ABP |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | influenza virus NS1A binding protein;IVNS1ABP;DKFZp686K06216;FLARA3;FLJ10069;FLJ10411;FLJ10962;FLJ35593;FLJ36593;HSPC068;KIAA0850;ND1;NS-1;NS1-BP;NS1BP ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IVNS1ABP |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |