



JSRP1 Antibody

Product Code	CSB-PA011969GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MG2
Immunogen	Human JSRP1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	junctional sarcoplasmic reticulum protein 1;JSRP1;FLJ32416;JP-45 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	JSRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.