



KANK2 Antibody

Product Code	CSB-PA011983GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63ZY3
Immunogen	Human KANK2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	KN motif and ankyrin repeat domains 2;KANK2;ANKRD25;DKFZp434N161;FLJ20004;KIAA1518;MGC119707;MXRA 3;SIP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KANK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.