



KBTBD2 Antibody

Product Code	CSB-PA011997GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY47
Immunogen	Human KBTBD2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	kelch repeat and BTB (POZ) domain containing 2;KBTBD2;BKLHD1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KBTBD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.