



KDM3B Antibody

Product Code	CSB-PA012136GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7LBC6
Immunogen	Human KDM3B
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	lysine (K)-specific demethylase 3B; KDM3B; 5qNCA; C5orf7; JMJD1B; KIAA1082 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KDM3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.