



KIF1B Antibody

Product Code	CSB-PA012319GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60333
Immunogen	Human KIF1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	kinesin family member 1B;KIF1B;CMT2;CMT2A;CMT2A1;FLJ23699;HMSNII;KIAA0591;KIAA1448;KLP ;MGC134844 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KIF1B
Target Details	This gene encodes a motor protein that transports mitochondria and synaptic vesicle precursors. Mutations in this gene cause Charcot-Marie-Tooth disease, type 2A1.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.