



KLHL29 Antibody

Product Code	CSB-PA012430GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CT2
Immunogen	Human KLHL29
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	kelch-like 29 (Drosophila);KLHL29;KBTBD9;KIAA1921 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KLHL29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.