



KMO Antibody

Product Code	CSB-PA012475GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15229
Immunogen	Human KMO
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	kynurenine 3-monooxygenase (kynurenine 3-hydroxylase);KMO;dJ317G22.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KMO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.