



L3MBTL3 Antibody

Product Code	CSB-PA012712GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JM7
Immunogen	Human L3MBTL3
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	I(3)mbt-like 3 (Drosophila); L3MBTL3; MBT-1; MBT1; RP11-73O6.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	L3MBTL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.