



LAMP2 Antibody

Product Code	CSB-PA012740GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13473
Immunogen	Human LAMP2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	lysosomal-associated membrane protein 2; LAMP2; CD107b; LAMPB; LGP110 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LAMP2
Target Details	This protein is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.