



# Recombinant Mouse Transcription factor SOX-1 (Sox1)

<b>Product Code</b>	CSB-EP022417MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53783
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p> MYSMMMETDL HSPGGAQAPT NLSGPAGAGG GGGGGGGGGG  GGGTKANQDR VKRPMNAFMV WSRGQRRKMA QENPKMHNSE  ISKRLGAEWK VMSEAEKRPF IDEAKRLRAL HMKEHPDYKY RPRRKTCTLL  KKDKYSLAGG LLAAGAGGGG AAVAMGVGVG VGAAAVGQRL  ESPGGAAGGG YAHVNGWANG AYPGSVAAAA AAAAMMQEAQ  LAYGQHPGAG GAHPHAHPAH PHPHHPHAHP HNPQPMHRYD  MGALQYSPIS NSQGYMSASP SGYGGIPYGA AAAAAAAGG AHQNSAVAAA  AAAAAASSGA LGALGSLVKS EPSGSPAPA HSRAPCPGDL REMISMYLPA  GEGGDPAAAA AAAAQ SRLHS LPQHYQGAGA GVNGTVPLTH I </p>
<b>Source</b>	E.coli
<b>Target Names</b>	Sox1
<b>Protein Names</b>	Recommended name: Transcription factor SOX-1
<b>Expression Region</b>	1-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. In mice, a similar protein regulates the gamma-crystallin genes and is essential for lens development.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.