



# LDHC Antibody

<b>Product Code</b>	CSB-PA012844GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P07864
<b>Immunogen</b>	Human LDHC
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,WB,IF
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	lactate dehydrogenase C;LDHC;CT32;LDH3;LDHX;MGC111073 ;
<b>Product Type</b>	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	LDHC
<b>Target Details</b>	Lactate dehydrogenase C catalyzes the conversion of L-lactate and NAD to pyruvate and NADH in the final step of anaerobic glycolysis. LDHC is testis-specific and belongs to the lactate dehydrogenase family. Two transcript variants have been detected which differ in the 5 untranslated region.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.