



# Recombinant Human Calcipressin-1 (RCAN1)

<b>Product Code</b>	CSB-BP019500HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53805
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEEVDLQDLPSATIACHLDPRVFDGLCRAKFESLFRTYDKDITFQYFKSFKRV RINFSPFSAADARLQLHKTEFLGKEMKLYFAQTLHIGSSHLAPPNDKQFLIS PPASPPVGWKQVEDATPVINYDLLYAISKLGPGKEYELHAATDTTPSVVVHVC ESDQEKEEEEEEMERMRRPKPKIIQTRRPEYTPIHLS
<b>Source</b>	Baculovirus
<b>Target Names</b>	RCAN1
<b>Protein Names</b>	Recommended name: Calcipressin-1 Alternative name(s): Adapt78 Down syndrome critical region protein 1 Myocyte-enriched calcineurin-interacting protein 1 Short name= MCIP1 Regulator of calcineurin 1
<b>Expression Region</b>	1-197
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length of isoform 2
<b>Target Details</b>	This protein interacts with calcineurin A and inhibits calcineurin-dependent signaling pathways, possibly affecting central nervous system development. This gene is located in the minimal candidate region for the Down syndrome phenotype, and is overexpressed in the brain of Down syndrome fetuses. Chronic overexpression of this gene may lead to neurofibrillary tangles such as those associated with Alzheimer disease. Three transcript variants encoding three different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.